

# **A 1 a - Powder Moisture Accurate Standard Method**

## **GEA NIRO® Method No. A 1a**

Revised: January 2024

### **1. Definition**

The moisture content of a powder is the loss in weight (%) after oven drying at 102°C until constant weight is obtained.

### **2. Scope**

This is an accurate standard method which may be used for milk powder and all other powdered dairy products which do not contain crystallized lactose ( $\alpha$ -lactose-monohydrate).

### **3. Principle**

The sample is dried by oven drying to constant weight at  $102^{\circ}\text{C} \pm 2^{\circ}\text{C}$  for 2 hours. The oven drying is repeated until the two successive weighings do not differ more than 0.5 mg.

### **4. Apparatus**

- 4.1 Drying oven, with thermostat and without forced air circulation.
- 4.2 Analytical balance, sensibility  $\pm 0.1$  mg.
- 4.3 Desiccator with colour-indicating desiccant (e.g. silica gel).
- 4.4 Weighing dishes with lid.

### **5. Reagents**

None

### **6. Procedure**

- 6.1 Dry weighing dish with open lid in the oven and cool it in desiccator.
- 6.2 Weigh the empty dish (a), add approx. 3 g of powder and weigh again (b)
- 6.3 Place the loaded dish with open lid in the oven at  $102^{\circ}\text{C} \pm 2^{\circ}\text{C}$  for 2 hours.
- 6.4 Cool closed dish to room temperature in desiccator, and weigh (c).
- 6.5 Continue drying the loaded dish with open lid in the oven at  $102^{\circ}\text{C} \pm 2^{\circ}\text{C}$  for 1 hour.
- 6.6 Repeat the cooling 6.4 and weigh again (c).

6.7 Repeat 6.5 until weight (c) is constant (i.e. until two successive weighings differ less than 0.5 mg)

## 7. Calculation

$$\text{Moisture} = \frac{b-c}{b-a} \times 100\%$$

a = weight of empty dish

b = weight of dish + powder

c = weight of dish + dried powder

## 8. Reproducibility

± 0.1 %

## 9. Remarks

A sample for moisture determination has to be handled carefully in order to avoid evaporation or prevent adsorption.

## 10. Literature

- GEA Niro Research Laboratory
- IDF Standard № 26:2004 / ISO Standard № 5537:2004
- De Knecht, R.J. and Brink, H.v.d.: Improvement of the drying oven method for the Determination of the Moisture Content of Milk Powder. Int. Dairy Journal, 8, 1998, pp. 733-738.